

Hydraulic Bed



Achievements & Awards :



TUV SUD - ISO 9001



Corporate Branding Award



Platinum



Enterprise



SMB - Maybank Rising Star



SMB Emerging Business



TUV Product Service

Hospital Hydraulic Bed

BA 4000 series






The BA 4000 series was developed with the healthcare personnel and patients in mind. It provides a safer and smarter alternative to regular hospital beds. Its advanced hydraulic and gas spring system allows care givers to operate the bed with ease, hence eliminating occupational hazards and enhancing the overall quality of patient care.

The BA 4000 series... a smarter choice...

- ① Hydraulic Bed (BA 4120-T)
- ② Bedside Rack (WP 1100)
- ③ Overbed Table (WP 2200)



... Beds utilize only the highest quality
"enabling technology" components.
... on the relevant standards and



VDA




Hospital Hydraulic Bed

BA 4000 series



Features

- Hydraulic foot pedal controlled Hi-Lo
- Gas spring controlled backrest and legrest
- Gas Spring trendelenburg/reverse trendelenburg tilting
- Detachable ABS head and foot boards
- 150mm central locking wheels




BA 4120 - T

- Material : Steel Framework
- Finishing : Scratch and chip resistant epoxy coating



Features

- Hydraulic foot pedal controlled Hi-Lo
- Gas spring controlled backrest and retractable crank legrest
- Gas Spring trendelenburg/reverse trendelenburg tilting
- Detachable ABS head and foot boards
- 150mm central locking wheels



BA 4021 - T

- Material : Steel Framework
- Finishing : Scratch and chip resistant epoxy coating

Bed Positions

- Backrest and Legrest (0° - 80°) (0° - 40°)

- Hi-Lo Position (460 mm - 820 mm)

- Trendelenburg and Reverse Trendelenburg (0° - 12°) (0° - 6°)



Hospital Hydraulic Bed Accessories

Product Features



Bed Stripper

For Linen or personal items placement.



Bed Tilting Levers

To achieve trendelenburg and reverse trendelenburg position.



Gas Spring Lever

This system allows fowler positions to be adjusted with greater ease.



Castors

150 mm diameter high quality German PU castors.

Optional



BX 200

Collapsible Aluminium Siderail

- Removable and secured by screw knob.
- Easily raised or lowered.



BX 120

Telescopic Siderail (Chrome)

- Chrome plated steel construction
- Extendable length to accommodate various bed sizes (auto locking).
- Also available in Stainless Steel (BX 121)



Central Locking Castors

150 mm diameter Tente (Germany) twin castors with unrivalled durability and mobility.

Hospital Hydraulic Bed

BA 4000 Series Beds

	Head/Foot Board		Fowler Position		Fowler Mechanism			Special Features	
	Wood panel with steel frame	Poly urethane	Single	Double	Head	Leg		Trendelenburg	X-Ray Translucent Base
					Gas Spring	Crank	Gas Spring		
BA 4010	●		●		●				
BA 4020	●			●	●		●		
BA 4021	●			●	●	●			
BA 4020-T	●			●	●		●	●	
BA 4021-T	●			●	●	●		●	
BA 4020-TX	●			●	●		●	●	●
BA 4110		●	●		●				
BA 4120		●		●	●		●		
BA 4121		●		●	●	●			
BA 4120-T		●		●	●		●	●	
BA 4121-T		●		●	●	●		●	
BA 4120-TX		●		●	●		●	●	●

**All rights reserved by LKL Advance Metaltech Sdn. Bhd.
Due to technological and design innovations, specifications of products may be subjected to change without prior notice.*



LKL Advance Metaltech Sdn.Bhd. (278577-X)
 29-31, 33-35, Jalan Indah 2/16, Taman Universiti Indah,
 43300 Seri Kembangan, Selangor Darul Ehsan, Malaysia.
 Tel: 603-8942 5230, 8948 2990, 8943 2050, 8943 5727
 Fax: 603-8948 7904, 8943 8487
 e-mail: info@lklbeds.com
 www.lklbeds.com