

Hospital Electrical Bed



The *Excellent Guide*

Towards Patient Care



www.lklbeds.com



Hospital Electrical Bed - BA7000 Series

BA7001 Ultra Low Bed



:: Features

- + Ultra low with the minimum height at 240mm (9.5")
- + Four section bed
- + All parts are made by tooling to ensure uniform quality
- + Finishing with epoxy polyester coated & baked finish resistant to chipping & scratching
- + Electric motorize bed with smooth and quite operation by 3 motors.
- + Detachable head and foot boards
- + Mobile on 100mm diameter castors with 2 brakes diagonally fitted
- + Provide with drainage bag holder at both side of the bed
- + Hilo, backrest & knee controlled by hand set
- + Provision for IV pole
- + With rotating corner bumper

:: Electrical System

- + Linak DC-actuator system (TUV & UL approved & EC60601-1)
- + 24 DC Output
- + Voltage: 220-240; 50Hz - 60Hz

:: Technical Specification

- + Mattress Dim(mm): 1980(L) X 864(W) x 102(H)

:: Optional

- + Monkey/Lifting Pole (BX910)



Backrest: 0° - 75°



Kneerest: 0° - 45°



2190mm

1013mm



650mm

BA7002-HC Electrical Bed



:: Features


- + Four section bed
- + All parts are made by tooling to ensure uniform quality
- + Steel base bed-tops & frame are made with epoxy powder coated & baked finish resistant to chipping and scratching
- + Electric motorize homeware bed with smooth and quiet operation
- + Detachable PE head / foot boards
- + Four high-quality 125mm swivel with 2 brakes
- + Hilo, backrest & knee controlled by hand set

:: Technical Specification

- + Working capacity: 180kg

:: Electrical System:

- + LINAK DC-actuator system (TUV & UL approved)
- + Voltage: 220 - 240; 50Hz - 60Hz



BA7003 Electrical Bed



:: Features


- + Four section bed
- + All parts are made by tooling to ensure uniform quality
- + Finishing with epoxy polyester coated & baked finish resistant to chipping & scratching
- + Electric motorize homecare bed with smooth and quiet operation
- + Detachable head and foot boards
- + Mobile on 150mm diameter castors with centre locking sytem with steer, free and brake.
- + Provide with drainage bag holder at both side of the bed
- + Provision for IV pole
- + Auto contour backrest and legrest adjust simultaneously
- + With rotating corner bumper

:: Electrical System

- + Linak DC-actuator system (TUV & UL approved & EC60601-1)
- + 24 DC Output
- + Voltage: 220 - 240; 50Hz - 60Hz
- + Built in back-up battery

:: Technical Specification

- + Working capacity: 180kg



Optional
for Model BA7003



Trendelenburg: 0° - 12°
Reverse Trendelenburg: 0° - 12°
(Model BA7003-T)

BA7733

Electrical Bed

With Manual Override



:: Features

- + Four section bed
- + All parts are made by tooling to ensure uniform quality.
- + Finishing with epoxy polyester coated & baked finish resistant to chipping & scratching.
- + With adjustable height (HiLo position), backrest, kneerest, by electrical & manual override by cranking lever at foot end of bed
- + Detachable head and foot boards
- + Mobile on 125mm diameter castors with 2 brakes diagonally fitted
- + With round bumpers at 4 corners

:: Technical Specification

- + Working capacity: 180kg

:: Electrical System

- + LINAK DC-actuator system (TUV & UL approved & EC60601-1)
- + 24DC out-put
- + Voltage: 220-240V; 50Hz-60Hz




Manual Override



Backrest, Kneesest & HiLo adjustable

Electrical

Auto contour not available for the Model BA7733



BA7770-T Basic ICU Bed



:: Features

- + Four section bed
- + Include manual trendelenburg and reverse trendelenburg
- + All parts are made by tooling to ensure uniform quality.
- + Steel base bed-tops & frame are made with epoxy powder coated & baked finish resistant to chipping and scratching.
- + Electric motorize bed with smooth and quite operation by 3 motors.
- + Detachable head and foot boards.
- + Mobile on 150mm diameter castors with centre locking sytem with steer, free and brake.
- + Hilo, backrest & kneerest controlled by hand set
- + Auto contour backrest and legrest adjust simultaneously
- + Provision for IV pole
- + Drainage bag holder

:: Electrical System

- + LINAK DC-actuator system
(TUV & UL approved & EC60601-1)
- + 24V out-put
- + Voltage: 220-240V; 50Hz-60Hz
- + Built in back-up battery

:: Technical Specification

- + Working capacity: 180kg

:: Optional Accessories

- + CPR Release Lever (Model BA7770-T + CPR)
- + With X-Ray Translucent base from head to umbilicus with cassette holder (Model BA 7770-TX)
- + Memory foam mattress (BX417)
- + IV pole with 2 hooks - chrome plated



Manual Trendelenburg & Reverse Trendelenburg



150mm swivel central locking castors
(Steer, Free & Brake)



Trendelenburg: 0° - 12°
Reverse Trendelenburg: 0° - 6°



Backrest: 0° - 75°




Kneerest: 0° - 45°

Inclusive Accessories

BX411

MediFlex Mattress (PU Foam)



- + Fire retardant
- + Anti-bacterial
- + Anti-static
- + Waterproof
- + Density: 32D
- + Dim(mm): 1905(L) x 864(W) x 102(H)

Handset

Electrical HiLo, Backrest, Kneerest & Auto contour adjustable



BX200

Collapsible side rail



Optional Accessories



Chrome/Stainless steel plated

Epoxy powder coated

WM223

IV pole height adjustable



278577-X

LKL Advance Metaltech Sdn. Bhd.
 No.3, Jalan BS 7/18,
 Kawasan Perindustrian Bukit Serdang,
 Seksyen 7, 43300 Seri Kembangan,
 Selangor Darul Ehsan.
 Malaysia.
 Tel: +603 8948 2990 (Hunting Line)
 Fax: +603 8948 7904, 8943 8487
 Email: info@lklbeds.com
 URL: www.lklbeds.com

Achievements & Awards :



MPC Productivity Award
2007



Woman Entrepreneur
2007



MMEBE Platinum
Brand Award



MMEBE Corporate Branding
Export Market (Gold)



5th Malaysia
100 Outstanding SMEs



50th Enterprise Award
Ranked No.6



SMB Emerging Business



SMI - Maybank Rising Star