

HOSPITAL BED W. M. F. DOUBLE FOWLER
MODEL BA 2020-P

- Features:**
- Steel base of 1.2mm thick with ventilation holes
 - With adjustable backrest & kneerest by means of winding mechanism at foot end of bed
 - With easy detachable PE blow head & leg boards.
 - Built-in bumper at 4 corners to prevent knocking against walls
 - Mobile on 4 x 100mm swivel grey rubber wheels diagonally fitted with brakes.
 - Provision for IV drip pole attachment

Finishing: - Hard epoxy powder coating cream colour.

Dimension : 2060mmL x 955mmW x 570mmH. (without mattress)

Mattress size : 1905mmL x 864mmW x102mmH

