

BED, FOWLER DOUBLE FIXED HEIGHT C/W MATTRESS
MODEL BA 2110-TGS

- Features:**
- With adjustable backrest by gas spring / gas lift.
 - With easy detachable polymer (PU) head & leg boards.
 - Built-in bumper at 4 corners to prevent knocking against walls
 - Mobile on 4 x 125mm swivel grey rubber wheels diagonally fitted with brakes.
 - Provision for IV drip pole attachment
 - Instant 2 way tilt from 0 degree to 12 degree by a fail-safe friction lock & 6° foot down tilt by a fail safe friction lock with every 3° position.

Finishing: - Finish of hard epoxy powder coating (beige color)

Dimension: - 2060mmL x 915mmW x 560mmH. (Without mattress)



HOSPITAL BED - COLLAPSIBLE SIDE RAILS

MODEL BX 200



HOSPITAL MATTRESS

MODEL BX 410



I.V. POLES

MODEL BX 300

