

HOSPITAL MECHANICAL FITTING – HILO (DOUBLE FOWLER)
MODEL BA 3020-PFLSD



Picture show the full length side rail for BA3020-PFLSR

Features:

- Steel base of 1.2mm thick with ventilation holes
- With adjustable Height (Hilo position), Backrest and Legrest by means of 3 cranks winding mechanisms at foot end of bed with retractable handle.
- The mechanism of crank is a solid metal rod with gear through the solid metal ball bearing by casting.
- With easy detachable PE blow head & leg boards.
- Built-in bumper at 4 corners to prevent knocking against walls
- Full Length Side Rail
- Mobile on 4 x 125mm PU castors diagonally fitted with brakes
- Provision for IV drip pole attachment

Finishing:

- Hard baked epoxy powder coating

Dimension: 2090mmL x 955mmW

Hilo Position: 380mm to 710mm (without mattress)

Mattress size: 1905mmL x 864mmWx102mmH

Optional:

- Central locking castor