

HOSPITAL HYDRAULIC BED
MODEL BA 4120-T

- Features:**
- Metal base with ventilation holes – mattress platform.
 - Fowlers position adjustable by gas spring.
 - Height adjustable by foot pump hydraulic system.
 - Instant 2 way tilt from 0 degree to 12 degree by a fail-safe friction lock & 6° foot down tilt by a fail safe friction lock with every 3° position.
 - With detachable Plastic PU head & foot boards.
 - Bumper guard at 4 corners to prevent knocking against wall.
 - Provision for IV drip pole attachment.
 - Bed stripper provided.
 - Four 150mm with central locking systems.

Finishing: Hard epoxy powder coating – Beige colour.

Hilo position from 460mmH to 820mmH (without mattress)

Dimension Size: 2060mmL x 925mmW x 460mmH

Mattress size: 1920mmL x 864mmW

