

HOSPITAL HYDRAULIC BED
MODEL BA 4120-TX

Features: Metal base with ventilation holes – mattress platform.
Fowlers position adjustable by either crank mechanism / gas spring.
Height adjustable by foot pump hydraulic system.
Instant 2 way tilt from 0 degree to 12 degree by a fail-safe friction lock & 6° foot down tilt by a fail safe friction lock with every 3° position.
With detachable Plastic PU head & foot boards.
With X-Ray Translucent base from head to umbilicus with cassette holder.
Bumper guard at 4 corners to prevent knocking against wall.
Provision for IV drip pole attachment.
Bed stripper provided.
Four 125mm swivel ball bearing rubber wheels with central locking.

Finishing: Hard epoxy powder coating – Beige colour.

Hilo position from 460mm to 820mm (without mattress)

Mattress size: 192cmL x 84cmW

