

LKL – HOSPITAL ELECTRIC ICU/CCU BED

MODEL: BA6000

FEATURES

- All parts are made by tooling to ensure uniform quality
- Mattress base tops made of ABS (high impact plastic) easy to clean.
- Bed frame are made of steel with epoxy powder coated & baked finish resistant to chipping and scratching
- Easy detachable & removable and can clean easily Head & Foot boards (Material: PE blow in one piece).
- Four high-quality 125mmdia. Central locking system.
- Hi-Lo, Backrest, Legrest and Trendelenburg and Reverse Trendelenburg by electric motorize by patient control and nurse control.
- 4 pcs foldable PE blow guard rails , 2 pcs at each side of the bed with smooth and easy to clean surface.
- Able to get chair position.
- Provision of IV pole attachment at 4 corners of bed.
- Provision for monkey pole attachment
- Complete with PU foam Mattress with chemsafe cover (fire retardant, water proof, anti-bacteria and anti-static)

ELECTRIC SYSTEM

- LINAK DC-actuator system (TUV & UL approved & EC60601-1)
- 24V battery back-up function-optional
- Voltage: 220-240V; 50Hz-60Hz
- Back-up battery

FUNCTIONS

- Backrest & knee lifting are 0° -70° & 0° - 35° respectively
- Bed height adjustable 290mm
- Trendelenburg & Reverse Trendelenburg (15°)

Overall size: 2200mmL x 927mmW x 450mm to 740mm H.

Mattress Size: 1980mmL x 915mmW x 100mmH

OPTIONAL ACCESSORIES:

- Stainless steel IV Drip Pole with 2 hook
- X-Ray radio translucent base from head to torso
- Drainage hook
- Oxygen cylinder holder
- Monkey Pole
- Nurse lock for the hand set located at foot end of the bed



LKL Advance Metaltech Sdn. Bhd. (278577-X)

Hospital Electrical ICU/CCU Bed

BA 6000

- * ABS Mattress Base
- * Nurse Lock
- * Chair Position
- * Fully Motorize System

(Backrest, Legrest, Hi-Lo, Trendelenburg & Reverse Trendelenburg)



	User Friendly "Handset"	
Adjustable Backrest & Legrest		Nurse Lock
	ABS Mattress Base	

