

LKL ELECTRICAL BED
MODEL BA 7770-T

Features:

- ***Include manual trendelenburg and reverse trendelenburg***
- All parts are made by tooling to ensure uniform quality.
- Steel base bed-tops & frame are made with epoxy powder coated & baked finish resistant to chipping and scratching.
- Electric motorise bed with smooth and quite operation by 3 motors.
- Detachable PE head / foot boards.
- Collapsible type guard rails (detachable), optional with telescopic side rail (chrome)
- Four high-quality 150mm swivel with central locking wheels

Electric System:

- LINAK DC-actuator system (TUV & UL approved & IEC60601-1)
- 24V out-put
- Voltage: 220-240V; 50Hz-60Hz
- Built in back-up battery

Functions:

- Hilo, backrest & knee controlled by hand set
- Backrest & kneerest are 0° - 75° & 0° – 45° respectively & adjustable simultaneously
- ***Manul Trendelenburg(12°) & Reverse Trendelenburg(6°) assisted by gas spring***
- Height adjustable: 460mm to 770mm
- Mattress size: 1905mmL x 865mmW x 102mmH
- Overall size: 2050mmL x 1030mmW

ACCESSORIES INCLUDED

- Collapsible side rail
- Mediflex Mattress with chemsafe cover (fire retardant, water proof, anti-bacteria and anti-static) and high density (32D) PU foam. (BX411).

Optional

- CPR Release Lever

