

HOSPITAL ICU/CCU ELECTRICAL BED ~ INTEGRATED WITH WEIGHING SCALE
MODEL BA9900-WA



FEATURES:

- All parts are made by tooling to ensure uniform quality
- Sheet metal mattress base for durability with X-ray translucent base for the head to umbilicus section
- X-ray cassette holder tray with adjustable film positioning on the backrest section
- Steel bed frame with epoxy polyester powder coated & baked finish is resistant to chipping and scratching
- PE Head & Foot boards with locking mechanism for easy attachment detachment
- Tente twin castors 150mm diameter central locking and steer system with pedals located at the foot end for convenient operation
- Split fold down PE blow safety side rail with convenient one hand operation for engaging and disengaging side rail
- Height, Backrest, Legrest and Trendelenburg and Reverse Trendelenburg by electric motorize
- Embedded dual sided controls on the head section side rails – full function on the outer side for caregivers and limited function on the inner side to allow patient to have own desired position
- Auto contour function to reduce the chances of patient sliding to the foot end of bed when – to achieve the most ergonomic and comfortable ‘cardiac chair position’, while reducing the need to move patients to alternative seating
- Caregiver control (nurse lock and nurse control) located at foot end of the bed
- Nurse lock to lock out the individual articulation function of the side rail
- CPR one press button – to return the bed to the lowest position for height, backrest and legrest
- Manual Quick CPR for lowering backrest release lever at both sides of the bed
- Digital weighing scale with segment LED display and control panel on the foot board
- Zero scale function
- 6 IV pole mounting and monkey pole mounting

- IV drip pole storage
- Drainage bag holder located on both sides of the bed
- PU bumper at 4 corners

Physical :

- Backrest : 0° – 70°
- Knee rest : 0° – 35°
- Height : 480mm to 780mm (to bed pan level)
- Trendelenburg & Reverse Trendelenburg : 15°
- Overall size: 2200mmL x 1050mmW
- Sleeping surface : 1980mmL x 915mmW
- Safe load capacity : 230kg

Electrical system :

- Input Voltage: 220-240V; 50Hz-60Hz
- LINAK DC-actuator system (TUV & UL approved & EC60601-1)
- 24VDC with built in back-up battery
- Safety rating
 controllers : IP54
 actuators : IP54
 nurse control : IP66

Accessories included :

- PatientMate Technology Mattress with Memory Foam - allowing patient's body contouring and reducing pressure, cover of fire retardant, anti bacterial and fluid resistant, thickness 100mm
- Stainless Steel Telescopic IV Drip Pole with 2 hooks

Optional Accessories:

- Monkey Pole
- Oxygen Cylinder Holder
- Portable Monitor / Defibrillator Tray
- Orthopedic Frame