

HOSPITAL MATTRESS
MODEL BX 411

Mattress:

- 100% waterproof mattress cover
- Extra comfort: unique construction provides excellent stretch and recovery to prevent pudding
- sanitised : anti-bacteria (resists harmful bacteria)
- Urine resistant
- Stain retardant
- Fire retardant
- Anti-static
- Polymer type: Homopolymer polyvinyl chloride
- Colour: beige
- Thickness: 0.35 – 0.38mm
- Weight range: 380 – 414 gsm
- Hardness: 60phr
- Water proof property: applicable

Size: 75”L x 34”W x 4”H

Mattress foam:

D32 PU foam

