

HOSPITAL OT - OBSTETRIC BED
MODEL DB 1120

Made of steel frame constructed of 3 section s/steel base top

Tilting positions by a manual crank gear mechanism at side of table.

Head Section: Manually adjustable backrest by ratchet 4 positions.
PVC mattress of 75mm thick (2 hinge section) & a pair of lithotomy stirrups and straps. Detachable headboard and side rails. Inclusive of an IV drip pole.

Foot Section: Slides smoothly on steel bearings under the head section and locks automatically. PVC padded mattress 100mm thick - removable.

Castors: Four 125mm swivel ball bearing rubber wheels with brakes diagonally fitted.

Finishing: Hard epoxy powder coating - cream colour.

Dimension: 1180mm – 1830mmL (Extendable)x 760mmW x 750mmH.

** Available with a s/steel solution bowl. A removable step provided.

