

Fixed Height Delivery Bed



DB 1200

Delivery Bed


- + Steel framework with scratch and chip resistant epoxy powder coating
- + Mobile on 4 x 125mm dia. swivel castors with 2 brakes diagonally fitted
- + Removable siderails
- + Stainless steel I.V pole
- + Ratchet adjustable backrest
- + Stainless steel bowl
- + Removable foot step
- + Gas spring tilting
- + Hand grip
- + 2 lithotomy stirrups & straps
- + Mattress with PVC cover
- + Removable head board



Backrest
0° - 65°



180 kg
UDL



Trendelenburg : 0° - 12°
Rev. Trendelenburg : 0° - 12°



DB 1120

Delivery Bed

- + Steel framework with scratch and chip resistant epoxy powder coating
- + Mobile on 4 x 125mm dia. swivel castors with 2 brakes diagonally fitted
- + Removable siderails
- + Stainless steel I.V pole
- + Ratchet adjustable backrest
- + Stainless steel bowl
- + Removable head board
- + Tilting by crank
- + 2 lithotomy stirrups
- + Mattress with PVC cover
- + Removable foot step



Backrest
0° - 65°



180 kg
UDL

