

HOSPITAL OT - BIRTH CARE BED
MODEL DB 2000

Features:

- The multi-purpose Birth Care System Bed can be converted from a practical labour bed to a delivery platform and back to a comfortable recovery bed.
- It comes with an extensive set of practical accessories such as adjustable s/steel side rails, transfusion pole, fluid collection bowl, hand grips for bearing down, foot support, adjustable & comfortable leg rest for lithotomy position, padded 2-hinge mattress and a removable step for patient access.
- Instant head down tilting & raising backrest by Gas Spring Height adjustable by hydraulic system foot pedal at both side of bed.
- Special moulded lower platform completely removable for easy cleaning.
- Retractable leg section. Removable headboard for access by anesthetists.

Castors: available with central locking system with directional wheel by 2 foot operated pedals fitted at foot end.

Finishing: Hard epoxy powder coating - cream colour.

Dimension: 1990mmL x 810mm
Hilo position from 510mm to 880mm approx.

