

HOSPITAL OT – ELECTRIC BIRTH CARE BED
MODEL DB 3000

Features: The multi-purpose BirthCare System Bed can be converted from a practical labour bed to a delivery platform and back to a comfortable recovery bed. It comes with an extensive set of practical accessories such as adjustable s/steel side rails, transfusion pole, fluid collection bowl, hand grips for bearing down, adjustable & comfortable leg rest for lithotomy position, padded 2-hinge mattress and a removable step for patient access.

Retractable leg section.

Removable headboard for access by anaesthetist

ELECTRICAL system Height adjustment and Backrest adjustment

Instant head down tilting assisted by Gas Spring system.

Special moulded lower platform completely removable for easy cleaning.

Castors: With central locking system with directional wheel by foot operated pedals fitted at head end.

Finishing: Hard epoxy powder coating

Size: 1990mmL x 810mm
Hilo position from 510mm to 880mm approx.

