

HOSPITAL BED W. M. F. CHILD COT
MODEL DC 2010

Steel base of 1.2mm thick with ventilation holes.

Tubular steel pole head & leg frame

Safety side rails of filling not exceeding 4" spacing with rubber bushes along 2 chrome rods fitted on leg frame eliminating friction noises.

Engaged automatically when raised.

Single fowler back rest adjustable by ratchet.

Castors : Four 75mm swivel castors with brakes diagonally fitted.

Finishing : Hard epoxy powder coating cream colour.

Dimension : 1400mmL x 760mmW x 610mmH (without mattress)

Mattress size: 127cmL x 76cmW

