

Examination Table



Achievements & Awards :





Examination Table



EC 2001

Examination Table



- Steel construction
- Scratch and chip resistant epoxy coating
- PVC cushion
- Ratchet back rest
- 2 cabinets
- Dim(mm) :



EC 3000

Examination Table - Adjustable Height

- Made of tubular steel leg frame with sectional base (2 section)
- Finishing of hard epoxy powder coating BEIGE colour
- Adjustable backrest (0° to 90°) - assisted by gas spring with a control lever on both side of table
- Height adjustable by hydraulic system foot pedal at both side of table
- Cushion of PVC padded foam mattress cover 75mm thick
- Inclusive of a paper roll holder at head end
- Mobile on 4 x 125mm swivel castors with 2 brakes
- Dim(mm) :



* optional: paper roll, adjustable legrest, Telescopic Side Rail, Trendelenburg/reverse trendelenburg position



EC 4001

Examination Massage Table

- Tubular mild steel frame c/w heavy base stand
- Armrest & Top mattress grey PVC covered 50mm thick rubber foam
- Adjustable lowered down ratchet system backrest c/w centre hole -15°
- Finishing of hard epoxy powder coating GREY texture
- Dim(mm): 1778 (L) x 560 (W) x 585 (H)


Examination Table



EC 1000

Examination Table


- Steel construction
- Scratch and chip resistant epoxy coating
- PVC cushion
- Ratchet back rest
- Dim(mm) :



EC 1001

Examination Table


- Steel construction
- Scratch and chip resistant epoxy coating
- PVC cushion
- Ratchet back rest
- Foldable foot step
- Dim(mm) :



EC 2000

Examination Table

- Steel construction
- Scratch and chip resistant epoxy coating
- PVC cushion
- Ratchet back rest
- 1 cabinet and 2 drawers
- Dim(mm) :




Examination Table



EC 5001

Examination Table For Gynae

- Steel construction ● S. steel bowl
- Scratch and chip resistant epoxy coating ● PVC cushion
- 2 Lithotomy stirrups and straps ● Ratchet back rest and leg rest
- Dim(mm) :




*available in EC 5000 (without S. steel bowl)



EC 6000

Examination Table For Gynae


- Steel construction ● S. steel bowl ● 1 cabinet and 3 drawers
- Scratch and chip resistant epoxy coating ● PVC cushion
- 2 Lithotomy stirrups and straps ● Ratchet back rest and leg rest
- Dim(mm) :



EC 7000

Examination Table For Gynae - Adjustable Height

- Made of tubular steel leg frame with sectional base (3 section)
- Finishing of hard epoxy powder coating BEIGE colour
- Manual adjustable backrest by ratchet and cut-away/lowered-down footrest
- Height adjustable by hydraulic system foot pedal at both side of bed
- Cushion of PVC padded foam mattress cover 75mm thick
- Inclusive of a pair of adjustable legrest, a stainless steel tray and paper roll holder at head end
- Mobile on 4 x 125mm swivel castors with 2 brakes
- Dim(mm) :



*optional: paper roll

**All rights reserved by LKL Advance Metaltech Sdn. Bhd.*

Due to technological and design innovations, specifications of products may be subjected to change without prior notice.



LKL Advance Metaltech Sdn.Bhd. (278577-X)

29-31, 33-35, Jalan Indah 2/16, Taman Universiti Indah,

43300 Seri Kembangan, Selangor Darul Ehsan, Malaysia.

Tel: 603-8942 5230, 8948 2990, 8943 2050, 8943 5727

Fax: 603-8948 7904, 8943 8487

e-mail: info@lklbeds.com

www.lklbeds.com