

# Healthcare Electrical Bed




# Healthcare Electrical Bed BA 7002-HC



## Functions


- Hi-lo, Backrest and Kneerest controlled by hand set.



- Backrest & knee lifting are  $0^{\circ} - 80^{\circ}$  &  $0^{\circ} - 45^{\circ}$  respectively & adjustable simultaneously.




- Height adjustable:  
Lowest position: 360mm  
Highest position: 670mm



- Minimum height: 360mm




- Mattress size: 1905mmL x 914mmW x 102mmH




## Overall size

2110mmL x 970mmW



360mm to 670mmH






## Features

- \* All parts are made by tooling to ensure uniform quality
- \* Steel base bed-tops & frame are baked epoxy-polyester finishing resistant to chipping and scratching
- \* Electric Motorize patient bed with smooth and quiet operation
- \* Detachable head / foot board with PU blow (Plastic)
- \* Four high-quality 125mm PU castors with 2 brakes.

## Electric System

- \* LINAK DC-actuator system (TUV & UL approved)



## Inclusive Accessories

- \* Extruded aluminum collapsible type guard rails (detachable)

- \* Mediflex mattress: Fire retardant, water-proof PU foam Mattress



LKL Beds utilize only the highest quality "enabling technology" components. Among the relevant standards and certification bodies are :



- UL 2601-1



- DIN EN ISO 14001



## Healthcare Electrical Bed

\*All rights reserved by LKL Advance Metaltech Sdn. Bhd.  
Due to technological and design innovations, specifications of products may be subjected to change without prior notice.



LKL Advance Metaltech Sdn. Bhd.  
No. 29-31, 33-35, Jalan Indah 2/16,  
Taman Universiti Indah,  
43300 Seri Kembangan,

Tel : 603-8942 5230  
8948 2990  
8943 2050

