

HOSPITAL HILO A & E TROLLEY
MODEL PT2010

Features:

- Fixed Height Trolley. Tubular steel leg frame c/w mild steel base top.
- Detachable head panel. Manually adjustable backrest by ratchet 5 position.
- A utility tray & oxygen cylinder bracket provided below.
- Telescopic safety side rails.
- PVC padded mattress 50mm thick - 2 hinge section.
- Instant 2 way tilting by Gas Spring
- IV drip pole attachment & an IV double hook pole provided
- Protective moulded PVC bumper buffer at four corners to prevent knocking against walls

Castors: Four 200mm swivel ball bearing rubber wheels with top brakes diagonally fitted

Make: Hard epoxy powder coating – beige colour.
Tubular steel frame with manual raising backrest.

Dimension: 1830mmL x 610mmW (mattress)
2060mmL x 610mmW x 810mmH

