

## **HOSPITAL OT – ELECTRIC BIRTH CARE BED** **MODEL DB 4000**

**Features:** The multi-purpose BirthCare System Bed can be converted from a practical labour bed to a delivery platform and back to a comfortable recovery bed. It comes with an extensive set of practical accessories, transfusion pole, fluid collection bowl, hand grips for bearing down, adjustable & comfortable leg rest for lithotomy position and a removable step for patient access.  
3 section medisoft mattress (Blue Colour)  
PE side rails.  
Retractable leg section.  
Removable PE headboard for access by anaesthetist

**ELECTRICAL system Height, Backrest and trendelenburg adjustment**  
Special moulded lower platform completely removable for easy cleaning.

**Castors:** With 6” central locking system with directional wheel by foot operated pedals fitted at head end.

**Finishing:** Hard epoxy powder coating with “SAND BEIGE” Colour.

**Overall Size:** 2108mmL x 953mmW (full extend)  
Hilo position from 600mm to 990mm approx.



