



Arterial Venous Fistula Needle (AVF)

High Quality Performance with Extra Thin Wall Needle

Specifications

- Wing : Turnable Wing or Fixed Wing
- Cannula : Extra thin wall with / without back-Eye

- Tube Length : 300mm
150mm
- Sterilization : ETO / Gamma

Connector & Clamp

Color	Code
White	W
Red	R
Blue	B



Product Range

Gauge	Needle Length	Tube Length	Wing Color		Code	Back-Eye	Clamp	Wing		
			NIPRO	ISO				Turnable	Fixed	
14G	4/5" (20mm)	150mm or 300mm	Violet	Pale-Green	-				○	
	1" (25mm)				H	○			○	
	1-1/4" (32mm)				T		○			
15G	4/5" (20mm)	150mm or 300mm	Ivory	Blue*	Blue-Grey	C		○	○	
	1" (25mm)					HC	○	○		○
	1-1/4" (32mm)					HT	○		○	
16G	4/5" (20mm)	150mm or 300mm	Light-Green	White	TC		○	○		
	1" (25mm)				HTC	○	○	○		
	1-1/4" (32mm)									
17G	4/5" (20mm)	150mm or 300mm	Orange	Red-Violet						
	1" (25mm)									
	1-1/4" (32mm)									

*Single needle available with needle length 1" (25mm)

Package : 500 Pcs / Carton for 300mm Tube Length
1000 Pcs / Carton for 150mm Tube Length

* for USA market